

# Ritter Stools

## Comparison Chart

| Features   | Ritter Value Series Stools  |              |              |              | Ritter Classic Series Stools   |              |              |              |             |              |              |
|--|---|--------------|--------------|--------------|--|--------------|--------------|--------------|-------------|--------------|--------------|
|  | 270   | 271          | 272          | 273          | 274  | 275          | 276          | 277          | 278         | 279          |              |
| Seat height  | 18" - 24.5"   | 18" - 24.5"  | 18" - 24.5"  | 18" - 24.5"  | 18" - 24.5"  | 18" - 24.5"  | 18" - 24.5"  | 18" - 24.5"  | 18" - 24.5" | 18" - 24.5"  | 19" - 25.5"  |
| Seat diameter                                      | 16.25"  | 16.25"       | 16.25"       | 16.25"       | 16.25"   | 16.25"       | 16.25"       | 16.25"       | 16.25"      | 16.25"       | 16.25"       |
| Seat cushion thickness                             | 3.25"   | 3.25"        | 3.25"        | 3.25"        | 3.25"  | 3.25"        | 3.25"        | 3.25"        | 3.25"       | 3.25"        | 3.25"        |
| Weight capacity                                    | 300 lb  | 225 lb       | 300 lb       | 225 lb       | 350 lb   | 225 lb       | 350 lb       | 225 lb       | 225 lb      | 225 lb       | 225 lb       |
| Composite black 5-leg caster base and hood         | •   | •            | •            | •            |  |              |              |              |             |              |              |
| Polished aluminum/black 5-leg caster base and hood |   |              |              |              | •  | •            | •            | •            | •           | •            | •            |
| Caster base diameter                               | 23"   | 23"          | 23"          | 23"          | 21.5"  | 21.5"        | 21.5"        | 21.5"        | 21.5"       | 21.5"        | 21.5"        |
| Backrest   |   | •            |              | •            |  | •            |              | •            |             |              |              |
| Backrest dimension                                 |   | 14"W x 7.5"H |              | 14"W x 7.5"H |  | 14"W x 7.5"H |              | 14"W x 7.5"H |             | 14"W x 7.5"H |              |
| Armrests   |   |              |              |              |  |              |              |              | Surgical    | Surgical     |              |
| Height adjustment mechanism                        | Screw Adjust  | Screw Adjust | Air Lift     | Air Lift     | Screw Adjust   | Screw Adjust | Air Lift     | Air Lift     | Air Lift    | Air Lift     | Air Lift     |
| Height control                                     | Manual  | Manual       | Single lever | Single lever | Manual   | Manual       | 8" Hand ring | 8" Hand ring | Foot lever  | Foot lever   | 8" Hand ring |
| Upholstery color options                           | 27  | 27           | 27           | 27           | 27   | 27           | 27           | 27           | 27          | 27           | 27           |
| Optional accessories available                     | <b>Ritter Value Series Stools</b><br>Locking casters, Auto-locking casters, Soft rubber casters, Glides, 6" lab extension, Foot ring. |              |              |              | <b>Ritter Classic Series Stools</b><br>Locking casters, Auto-locking casters, Soft rubber casters, Glides (not available on 278), 6" lab extension, Foot ring. |              |              |              |             |              |              |

opt = optional

For product inquiries, please contact the Midmark Medical Customer Experience Department at 1-800-MIDMARK, extension 108921.  
 For more information or a demonstration, contact: Midmark Corporation, Dayton, OH.  
 Call: 1-800-MIDMARK Fax: 1-800-365-8631 midmark.com

MKT-00004

