



### ADVANCED LIGHTING DESIGN

Midmark lights are engineered to deliver the true, reduced-shadow white light you need to properly match shades, identify details and diagnose tissue. And the automatic on/off function and upfront touchpad control system make the lights easy to use. Midmark lights are available in both LED and halogen options as well as various mounting configurations, including the Track Light Monitor, to meet your unique needs.

- Energy efficient 22-Watt power consumption
- 24" monitor available
- 5-year limited warranty



#### LED OPERATORY LIGHT

Offers energy efficiency and delivers true, reduced-shadow white light that allows you to properly identify shades, see details and evaluate tissues.



#### HALOGEN OPERATORY LIGHT

Delivers value and reliability using white incandescent light that helps to properly identify shades, see details and evaluate tissues.



#### TRACK LIGHT MONITOR

Combines excellent dental lighting with highly flexible monitor positioning, enhancing the patient experience and easing patient acceptance of procedures.



#### MIDMARK 255 LED PROCEDURE LIGHT

Whether you practice in an office environment or specialized procedure area, your work requires superior, white, shadow-free illumination.

# LIGHTING SPECIFICATIONS

Dental Lights	LED	Halogen
Operate Mode Color Temperature	4,500 Kelvin 4,750 Kelvin 5,000 Kelvin	5,700 Kelvin
Operate Mode Intensity, Foot Candles (Lux)	1,400 fc (15,000 lx) 1,950 fc (21,000 lx) 2,500 fc (27,000 lx)	1,400 fc (15,000 lx) 1,950 fc (21,000 lx) 2,500 fc (27,000 lx)
Composite Mode Color Temperature	4,000 Kelvin	4,000 Kelvin
Composite Mode Intensity, Foot Candles (Lux)	550 fc (5,900 lx)	2,100 fc (2,2000 lx)
Auto On/Off	✓	✓
Operate (Composite-Safe) Mode	✓	✓
Composite Mode	✓	✓
Touchpad Controls	✓	✓
Low-Wattage Power Consumption	✓	
LED Technology	✓	

# MOUNTING OPTIONS



Left/Right



Console



Cabinet



Wall



Ceiling



Track



Track Light Monitor



Universal