

- 10-year limited warranty on kink valves
- 5-year limited warranty



DESIGNED FROM THE INSIDE OUT TO COMPLEMENT THE WAY YOU WORK

Whether you need a compact, economic delivery system like the Asepsis 21®, the integrated design of the Procenter or the cutting-edge technology of the Elevance®, you can rely on our proven technology. At the heart of each Midmark dental delivery system is our patented Kink Valve® technology – designed for an unprecedented level of reliability and control with minimal maintenance. Midmark delivery systems are also available in multiple configurations with accessories and options to meet your unique needs.



ELEVANCE INSTRUMENT DELIVERY SYSTEM

The luxurious styling and comfort of Elevance extends to this sophisticated delivery system, giving you all the finger-tip controls, technology and styling you expect. The Integration Control Module® manages all of your components from one location, giving you maximum control with simple operation.



PROCENTER INSTRUMENT DELIVERY SYSTEM

Thoughtful design, attention to detail and quality construction make up the Procenter dental delivery system. Integrated controls and switches are designed to provide simple, yet refined, access for precision and flexibility. The Procenter is designed to compliment the way you work while delivering exceptional reliability, available in chair-, cabinet- and wall-mounted systems.



ASEPSIS 21 INSTRUMENT DELIVERY SYSTEM

The Asepsis 21 provides quality and reliability in a sleek, compact, value-priced package. Surfaces are smooth and easy to clean, while air and water controls are located beneath a removable top to help protect them from contamination.

DELIVERY UNIT SPECIFICATIONS

Dental Chairs	Elevance	Procenter	Asepsis 21
Left-Right Delivery Option	✓	✓	✓
Console Delivery Option	✓	✓	✓
12:00 Rear Delivery Option		✓	✓
Continental Unit Available	✓		✓
Assistant's Unit Available	✓	✓	✓
5-Year Limited Warranty	✓	✓	✓
Kink-Valve® Technology (with 10-year limited warranty)	✓	✓	✓
Chair Control Option	✓	✓	✓
Handpiece Integration	✓	✓	
3-Way Syringe Included	✓	✓	✓
Air-Lock Flex Arm		✓	
Programmable Presets	✓		
Touch-Screen Controls	✓		
Needle-Valve Air/Water Adjustments		✓	✓

AVAILABLE ACCESSORIES:

(Not all available on all units.)

Handpiece Illumination System

Tubing and Lamp Fiber
Optic Handpiece
Tubing for ISO-C Fiber
Optic Handpiece

Ultrasonic Scalers

Dentsply Cavitron®
Steri-mate® 360
Ultrasonic Scaler (G139)
SATELEC® Pure Newtron®
SLIM Ultrasonic Scaler
SATELEC® Pure Newtron®
SLIM B.LED Ultrasonic
Scaler

Electric Handpiece System

Bien Air® MX2 LED, Single
Motor System
Bien Air® MX2 LED, Dual
Motor System

Curing Light

SATELEC® Mini LED
Curing Light

Intraoral Cameras

SOPRO® Camera USB2
Docking Station
SOPRO® 617 Camera
SOPRO® 717 First Camera
SoproLIFE® Camera
SoproCARE® Camera

CONTINENTAL DELIVERY OPTIONS:



Elevance Continental Delivery

For over-the-patient instrument delivery that combines advanced design concepts with ergonomic balance, look no further than the Elevance Continental Delivery System. Experience the difference this sophisticated unit makes on comfort for you and your patients.



Asepsis 21 Continental Delivery

The reliability of the Asepsis 21 extends to an over-the-patient delivery for keeping all of your handpieces and air and water connections in the most convenient location to complement the way you work.

Cavitron and Steri-mate are registered trademarks owned by Dentsply International, Inc.
Bien Air is a registered trademark owned by Bien-Air Holding SA.
SATELEC and Newtron are registered trademarks owned by SATELEC Corporation.
SOPRO, SoproLIFE and SoproCARE are registered trademarks owned by SOPRO Corporation.

007-10047-00 Rev. C1 (12/18)
© 2018 Midmark Corporation, Maimisburg, Ohio USA