

Midmark 255  
LED Procedure Light



# Midmark 255 LED – superior procedure lighting at an affordable price.

Designed to give you the quality you deserve with features that can help you provide the best care to your patients and your bottom line.

## Superior Light Design

Significantly more light intensity than similar products at 7,500 footcandles (over 80,000 lux), an excellent color rendering index of 95 out of 100, and a color temperature of 4,400K.

## Counterbalanced

Easily position the light precisely where you want it, and it stays there – no drifting!

## Unsurpassed Shadow Reduction

16 LEDs reflecting evenly off of 1,040 reflective facets will reduce the disruption of shadows while you work, producing an even distribution of light for uninterrupted illumination.

## Dual On/Off and Brightness Control

Turn the light on/off and adjust light intensity from controls on each side of the light head.

## Reduced Energy Consumption

LEDs have 100-times longer life than most halogen lamps while consuming less energy, resulting in significant savings year after year. (See chart)


## Cleaning and Care

Sleek contours, a sealed light head and an autoclavable handle help to save maintenance time and increase product life.



## The LED Difference

Lifetime Operating Cost (Including Bulb Replacement Cost)



● LED Light ● Halogen Light



Find us on  
**Facebook**

facebook.com/MidmarkAnimalHealth