



Recessed units not shown

QUICKCLEAN™ ULTRASONIC CLEANER

POWERFUL CLEANING – CONSISTENT RESULTS

Providing powerful, effective cleaning and consistent results, the QuickClean Ultrasonic Cleaner is a piece of equipment your dental office can't be without. The QuickClean's cutting-edge technology can help you create a safer, more efficient work environment by decreasing worker exposure to contaminants and potential sharps injuries while reducing the time and effort needed for cleaning.

LEAP

POWERFUL TECHNOLOGY

The QuickClean unit cleans thoroughly and consistently throughout the entire bath using Frequency-Leap technology. This advancement in cleaning efficiency helps ensure all of your instruments are fully cleaned the first time, every time, no matter where they are placed in the basket.



SIMPLE TO OPERATE

The QuickClean is easy to use right out of the box, so you can get your staff up and running with minimal training time. Features like a digital LCD control panel that displays prompts, options for preset or customized settings, intelligent software that remembers the last cycle settings, and a fluid sensor that ensures the bath always has sufficient water levels help simplify your processes and increase productivity.



RANGE OF SIZES AND DESIGN OPTIONS

The QuickClean comes in three table top sizes (1.2, 3.3 and 6.6 gal.) and two recessed options (3.3 and 6.6 gal.), so you can choose the option that best fits your configuration and workflow. And the stainless steel finish and attractive design coordinate well with any décor.

150%
MORE EFFECTIVE
THAN SONICLEAN

CUSTOMIZED
CYCLE TIMES FROM
1-99
MINUTES

FULL-SIZE
SAFETY BASKET
INCLUDED

3-YEAR
LIMITED
WARRANTY




Specifications, accessories and service part information on reverse.

SPECIFICATIONS:

QUICKCLEAN™ MODEL NUMBERS:	TABLE TOP UNIT: QC1-01	TABLE TOP UNIT: QC3-01 RECESSED UNIT: QC3R-01	TABLE TOP UNIT: QC6-01 RECESSED UNIT: QC6R-01
Capacity:	1.2 gal/4.5 L	3.3 gal/12.5 L	6.6 gal/25L
Table Top Units			
Overall Length:	13.5" (34.5 cm)	14.7" (37.5 cm)	21.6" (55 cm)
Overall Width:	10.2" (26 cm)	16.5" (42 cm)	16.5" (42 cm)
Overall Height:	10.4" (26.5 cm)	10.4" (26.5 cm)	14.3" (36.5 cm)
Tank Size:	11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)	13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)	19.9" L x 11.8" W x 7.9" H (50.5 x 30 x 20 cm)
Recessed unit cut out size:	N/A	15.55" x 14.37" (39.5 x 36.5 cm)	22.48" x 14.37" (57.1 x 36.5 cm)
Table Top Unit Weight:	12.7 lbs (5.76 kg)	18.6 lbs (8.42 kg)	29.6 lbs (13.44 kg)
Table Top Shipping Weight:	14.6 lbs (6.64 kg)	20.9 lbs (9.50 kg)	32.6 lbs (14.78 kg)
Recessed Unit Weight:	N/A	22.3 lbs (10.10 kg)	32.9 lbs (14.92 kg)
Recessed Unit Shipping Weight:	N/A	24.6 lbs (11.18 kg)	35.8 lbs (16.26 kg)
Operating Frequency:	32-38 kHz	32-38 kHz	32-38 kHz
Temperature up to:	140°F / 60°C	140°F / 60°C	140°F / 60°C
Stainless Steel Safety Basket:	Included	Included	Included
Electrical Requirements	Input: 115V +/- 10%, 60 Hz, 110W, 1.0A Output: 100W	Input: 115V +/- 10%, 60 Hz, 210W, 1.8A Output: 200W	Input: 115V +/- 10%, 60 Hz, 410W, 3.6A Output: 400W
Standards	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B	UL 61010-1:2012; CAN/CSA C22.2 No. 61010-1:2012; EN 60601-1-2:2007; EN 61326:2013; FCC CFR 47 Part 15 part B
Limited Warranty	3 Years	3 Years	3 Years

ACCESSORIES & SERVICE PARTS:



1| BEAKER ACCESSORIES

- 9A612001 - QC1 – 2 Beaker Accessory
- 9A613001 - QC3/QC3R – 4 Beaker Accessory
- 9A614001 - QC6/QC6R – 6 Beaker Accessory

2| EXTRA SAFETY BASKETS

- 002-10007-00 – QC1 Safety Basket
- 002-10008-00 – QC3/QC3R Safety Basket
- 002-10009-00 – QC6/QC6R Safety Basket

3| CLEANING SOLUTIONS

- 9A296001 Midmark General Purpose Cleaner (32oz)
- 9A297001 Midmark Tartar and Stain Remover (32oz)
- 9A298001 Midmark Enzymatic Cleaner (32oz)

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark
dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2017 Midmark Corporation, Dayton, OH.
Products subject to improvement changes without notice
Litho in U.S.A. 007-10007-00 Rev. C1 (1/18)

