



**153672-003 Bien Air Curing Light (Duo Delivery)**  
(includes items 1 thru 17)

Item	Description	Qty.
1	• Curing Light Box Assembly (includes items 2 thru 17)	1
2	•• Curing Light Kit	1
3	•• Velcro (loop)	2
4	•• Velcro (hook)	2
5	•• Check Valve	1
6	•• Uni-Clamp (1/8" OD Tube)	2
7	•• Fitting (1/8" Barb Tee)	2
8	•• Barb Fitting (#10-32 x 1/16")	1
9	•• Gasket	1
10	•• Fitting (1/8" Barbs x 10-32 CRS)	1
11	•• Wire Connector	2
12	•• Clear Tubing (1/4" OD x 1/8" ID)	AR
13	•• Clear Tubing (1/8" OD)	AR
14	•• Black Wire (18 GA)	AR
15	•• Red Wire (18 GA)	AR
16	•• Blue Wire (18 GA)	AR
17	•• Brown Wire (18 GA)	AR

Serial Numbers / Dates: all