



153632-007 Dual Digital Electric Handpiece System, (ProCenter Unit)
(includes items 1 thru 18)

Item	Description	Qty.
1	• Screw	4
2	• Digital Handpiece Cover	1
3	• Electric Handpiece Bracket	1
4	• Screw	3
5	• Standoff	3
6	• Bien Air Digital Electric Handpiece Kit	1
7	• Electric Handpiece Bracket	1
8	• Screw	2
9	• Velcro (hook)	1
10	• Velcro (loop)	1
11	• Fitting, Barbed Union (1/8")	1
12	• Check Valve	1
13	• Tubing, Clear (1/8" OD)	AR
14	• Fitting, Barbed (1/8")	1
15	• Gasket	1
16	• Fitting, Barbed Nipple (1/4")	1
17	• Tubing, Clear (3/16")	AR
18	• Pressure Switch (NO)	3

DA165800i

Serial Numbers / Dates: V87315 Thru V352592