



153632-009 Electric Handpiece System (Duo Delivery)
(includes items 1 thru 15)

Item	Description	Qty.
1	• Bien Air Electric Handpiece Kit	1
2	• Screw	4
3	• Electric Handpiece Bracket	1
4	• Velcro (loop)	3
5	• Velcro (hook)	3
6	• Fitting, Barbed Union, (1/8")	4
7	• Check Valve	2
8	• Tubing, (1/8" OD Clear)	AR
9	• Barb Fitting (1/8")	2
10	• Gasket	2
11	• Barbed Nipple Fitting	2
12	• Connector	4
13	• Tubing, (1/8" Violet)	AR
14	• Pressure Switch (NO)	3
15	• Tubing, Clear	AR

DA166001i

Serial Numbers / Dates: V87315 Thru V352592