



153765-002 Satelec Scaler, 12:00 FSC
(includes items 1 thru 14)
153765-006 Satelec Scaler, 12:00 FSC
(includes items 1 thru 14)

Item	Description	Qty.
1	• Accessories Wiring Harness	1
2	• Clear Tubing (1/8" OD)	AR
3	• Barb Fitting (1/16")	1
4	• Clear Tubing (1/4" OD)	AR
5	• Barb Fitting (1/16" to 1/8")	1
6	• Uni-Clamp	5
7	• Tee Fitting (1/16")	1
8	• Check Valve	1
9	• Pressure Switch	1
10	• Wire Connector	7
11	• Satelec Scaler Kit (w/o Tips)	1
12	• Velcro Tape (loop)	1
13	• Velcro Tape (hook)	1
14	• Wire	AR

DA179301i

Serial Numbers / Dates: all