


**Note**  
 For complete breakdown of handpiece tubing, refer to: Midmark Delivery Systems Service & Parts Manual (004-0282-00).

040-0010-109 **2**

015-1480-00 **3**



**4** 177398

**5** 177118

**153630-008 Handpiece Illumination System - 120 V**  
 (includes items 1 thru 6)

Item	Description	Qty.
1	• 120V Transformer (includes Terminals)	1
2	• Screw (new-style units - w/o nuts)	2
3	• Power Pack	1
4	• Orange Wire	AR
5	• Green/Yellow Wire	AR

DA163201i

**Serial Numbers / Dates:** all