

**Note**  
Additional wires and tubing not shown are included with kit.


**153631-003 Amdent Scaler System (PoCenter)**  
(includes items 1 & 15)

Item	Description	Qty.
1	• Check Valve	1
2	• Barb Fitting (#10-32 x 1/16")	1
3	• Gasket	1
4	• Barbed Nipple Fitting	1
5	• Clear Tubing (1/8" OD)	1
6	• Barb Fitting (1/8" to 1/4")	1
7	• Screw	4
8	• Amdent Holder Bracket	1
9	• Screw	2
10	• Scaler Cover	1
11	• EMS Amdent Scaler Assembly	1
12	• Screw (#8-32 x 3/16")	1
13	• Water Control Valve Knob	1
14	• EMS Amdent Scaler Wiring Diagram (includes item 10) [not shown]	1

14 015-1605-00

DA163700i

Serial Numbers / Dates: all



DA2024i

<b>Serial Numbers / Dates:</b>	all
--------------------------------	-----