



153631-004 Amdent Scaler System
(includes items 1)

Item	Description	Qty.
1	• Duo Delivery Scaler Assembly (includes items 2 thru 18)	Ref
2	•• Amdent Scaler Components	1
3	•• Velcro (hook)	2
4	•• Velcro (loop)	2
5	•• Barb Fitting (1/16" to 1/8")	1
6	•• Check Valve	1
7	•• Uni-Clamp (1/8" OD Tube)	4
8	•• Fitting (1/16" Barb Tee)	1
9	•• Clear Tubing (1/8" OD)	AR
10	•• Clear Tubing (1/4" OD)	AR
11	•• Scaler Wiring	1
12	•• Black Wire to Handpiece Tube	Ref
13	•• White Wire to Handpiece Tube	Ref
14	•• 1/8" OD Clear to Handpiece	Ref
15	•• Tubing to Drs. Instrument Head	Ref
16	•• Green Wire to Ground	Ref
17	•• Brown Wires to OV & 24V	Ref
18	•• Control Knob and Label	Ref

DA1638001

Serial Numbers / Dates: all