


Note

The Red, Orange, and Yellow tubing are from Scaler.
 Green tubing is from Handpiece Supply Air.
 Brown & Yellow wires are from Scaler.



153764-001 Cavitron Scaler, (ProCenter Unit)
 (includes items 1 thru 8)

Item	Description	Qty.
1	• Kit, Cavitron 4' Scaler Assy.	1
2	• Splice	2
3	• Plug (for 1/8" OD Tubing)	1
4	• Tubing, Clear (1/4" OD x 1/8" ID)	AR
5	• Barb Fitting (1/8" to 1/4" Tubing)	1
6	• Uni-Clamp (1/8" OD Tubing)	1
7	• Barb Fitting (1/8" to 1/8" Tubing)	1
8	• 3 Way Barb Tee (1/8" Tubing)	1

DA168803i

Serial Numbers / Dates: all