



Note
 The Red, Orange, and Yellow tubing are from Scaler.
 Green tubing is from Handpiece Supply Air.
 Brown & Yellow wires are from Scaler.

DA169001i

153764-002 Cavitron Scaler, (ProCenter Unit)
 (includes items 1 thru 7)

Item	Description	Qty.
1	• Kit, Cavitron 4' Scaler assy.	1
2	• Plug (for 1/8" OD Tubing)	1
3	• Tubing, Clear (1/4" OD x 1/8" ID)	AR
4	• Barb Fitting (1/8" to 1/4" Tubing)	1
5	• Uni-Clamp (1/8" OD Tubing)	1
6	• Barb Fitting (1/8" to 1/8" Tubing)	1
7	• 3 Wat Barb Tee (1/8" Tubing)	1

Serial Numbers / Dates:	all
--------------------------------	-----