



DA179401i

153765-003 Satelec Scaler Duo Delivery Unit
 (includes items 1 thru 9)
153765-007 Satelec Scaler Duo Delivery Unit
 (includes items 1 thru 9)

Item	Description	Qty.
1	• Clear Tubing (1/4" OD)	AR
2	• Barb Fitting (1/16" to 1/8")	1
3	• Uni-Clamp	1
4	• Clear Tubing (1/8" OD)	AR
5	• Check Valve	1
6	• Fitting, Barb Tee (1/16")	1
7	• Pressure Switch	1
8	• Wire Connector	1
9	• Satelec Scaler Kit (w/o Tips)	1

Serial Numbers / Dates: all