


LR Unit Arm Assembly		
Item	Description	Qty.
1	LR Unit Arm Assembly (includes items 2 - 10)	1
2	• Needle Bearing	NSS
3	• Slotted Sleeve Bearing	NSS
4	• Water Bottle Cover	1
5	• Unit Arm Assembly	NSS
6	• Water Bottle Quick Disconnect Cap Assembly (includes item 7)	1
7	•• Quick Disconnect Cap (includes O-ring 014-1076-52)	1
8	• Screw (Phillips Head)	2
9	• Set Screw	1
10	• Flange Bearing	NSS
11	Bearing Assembly	Ref
12	Bearing Cap	Ref
13	Screw (Hex Head)	Ref
14	Toggle Switch and Water Bottle	Ref
15	City Water Selector Switch and Water Bottle (not shown)	Ref

NSS = Not Sold Separately

Page Applies To:	Mounted to:	Delivery System Serial Numbers:
Delivery Units (LR)	<i>Elevance Chair</i>	<i>all</i>


029-4365-00 1a
 029-4365-02 1b



Air Toggle Switch and Water Bottle		
Item	Description	Qty.
1	Toggle Switch / Water Bottle and In-line Regulator a) With 1 L Water Bottle (includes items 2 thru 8) ..b) With 2 L Water Bottle (includes items 2 thru 9)	1
2	• Air Toggle Switch	1
3	• Water Bottle Kit a) 1 Liter Water Bottle b) 2 Liter Water Bottle	1
4	• Plug	1
5	• Uniclamps	6
6	• Yellow Tubing	AR
7	• Blue Tubing	AR
8	• In-line Air Regulator Kit (includes uniclamps)	1
9	2L Water Bottle Pickup Tube (029-4365-02 only)	1
10	Quick Disconnect Cap Assy in LR Unit Arm	Ref

AR = As Required

Page Applies To: <i>Delivery Units (LR)</i>	Mounted to: <i>Elevance Chair</i>	Delivery System Serial Numbers: <i>all</i>
---	---	--



029-4365-01 1a
029-4365-03 1b

002-1709-00 2

7 117479

6 002-1660-00

5 177530

002-0777-00 3a

122830 3b

177523 4

10

9

DA4283i

City Water Selector Switch and Water Bottle		
Item	Description	Qty.
1	City Water Selector Toggle Switch / Water Bottle a) With 1 L Water Bottle (includes items 2 thru 8) b) With 2 L Water Bottle (includes items 2 thru 9)	1
2	• City Water Selector Switch	1
3	• Water Bottle Kit a) 1 Liter Water Bottle b) 2 Liter Water Bottle	1
4	• Blue Tubing	AR
5	• Yellow Tubing	AR
6	• In-line Air Regulator Kit (includes uniclamps)	1
7	• Uniclamps	6
8	• Water Selector Switch Label	1
9	2L Water Bottle Pickup Tube (029-4365-03 only)	1
10	Quick Disconnect Cap Assy in LR Unit Arm	Ref

AR = As Required

Page Applies To: <i>Delivery Units (LR)</i>	Mounted to: <i>Elevance Chair</i>	Delivery System Serial Numbers: <i>all</i>
---	---	--