



Console Mounted Assistant's Arm		
Item	Description	Qty.
1	Console Mounted Assistant's Arm Assembly (includes items 2 - 8)	1
2	• Radius Arm	1
3	• Dowel Pin (Left Hand Shown - RH in Other Hole)	1
4	• Washer (1.75 O.D. x .06 Thick)	1
5	• Washer	2
6	• Bearing (.502 I.D. x .927 O.D.)	1
7	• Adapter Pivot	1
8	• Shoulder Screw (3/8"-16 x 1.125)	1
9	Bolt (Socket Head Cap Screw)	3
10	Set Screw	3

Page Applies To:	Mounted to:	Delivery System Serial Numbers:
Assistant's Units	<i>Elevance Console</i>	<i>all</i>